## Pet Pathogen Nucleic Acid Detection Kit

### USER MANUAL

#### 【Packaging】

PCR amplification reagent 8 tests  
User manual 1 piece

#### 【Intended Use】

This kit uses Fluorescence-based PCR analysis to detect pathogens in pleural effusion, ascites, stool, eye/nasopharyngeal swab, and blood specimens from canine, feline or other animals. It is applicable to the detection, diagnosis and epidemiological investigation of pet pathogens.

#### 【Methodology】

Use nucleic acid extraction kits to extract nucleic acids (DNA/RNA) from specimens. Under the action of highly efficient reverse transcriptase, cDNA strands complementary to RNAs are synthesized during the 1-step reaction, serving as DNA templates. When Taq DNA polymerase functions, copies of specific target sequences are amplified after cycles of denaturation, annealing and extension. Fluorescent probes hybridize to the amplifying target DNAs. The 5'→3' exonuclease activity of Taq DNA polymerase then causes the reporter and quencher groups of probes separate from each other, generating specific fluorescence. The fluorescent PCR instrument detects fluorescent signals, and analyzes the results based on the Ct values and the amplification curves.

#### 【Storage and shelf life】

Storage temperature 4°C-30°C; shelf life: 24 months, please use it within the expiration date.

#### 【Operating environment】

Recommended operating temperature is 20-25 °C.

#### 【Specimens】

Pleural effusion, ascites, stool swabs, anal swabs, eye/nasopharyngeal swabs, EDTA anticoagulant blood, urine.

#### 【Sampling requirements】

**Liquid sample:** (including blood sample, abdominal effusion, pleural effusion or other liquid samples, etc.)  
Take 200 μL of liquid sample into preservation buffer, and mix well. Then, transfer 200 μL of the mixture into the LB well for lysis. For urine sample, directly take 200 μL and add to LB well for lysis.

#### Swab sample:

1. Fecal swab: take proper amount with a swab
2. Anal swab: moisten the swab with dilution buffer, and then wipe specimen
3. Eye/nasopharyngeal swab: swab under the eyelid and fully wipe to collect eye swab samples; swab oral and nasal discharge properly to collect nasopharyngeal swab samples
4. After swab sample collection, quickly snap the swab handle and place in the preservation buffer, and then shake it to fully dissolve the pathogen on the swab into the preservation buffer.

It is recommended that nucleic acid extraction and detection should be performed immediately after sample collection. If storage needed, samples can be stored at 4°C temporarily. Samples need to be refrigerated during transportation.

#### 【Preparation before testing】

1. Please read the instruction manual carefully before testing, be familiar with each step, and strictly follow the requirements when using the kit.
2. Wear disposable gloves and masks, operate under the condition meeting the requirement of PCR testing environment.

**Before PCR amplification, nucleic acid extraction must be completed by extraction reagent**

#### PCR amplification: (using real time PCR analyzer)

1. Use a pipette to draw 20μL of nucleic acid solution from the EB well of the cartridge and add it to the PCR testing tube. Flick the tube with your finger to fully dissolve and mix the lyophilized reagent. Shake all the solution in the tube to the bottom (or centrifuge for 5-10 second) and put it into the sample well of the PCR analyzer.

#### Note:

- a) Do not suck up the bubble in the solution when pipetting the nucleic acid solution
  - b) The lyophilized reagent must be completely dissolved and mixed, and all the reaction liquid in the tube must be shaken to the bottom. Otherwise, strong background signal interference might happen easily, which will affect the result analysis.
2. Select the corresponding program for the specific pet pathogens, set the sample and test information, and then press the operation icon to start the reaction.
  3. When the reaction is completed, take out the PCR product tube immediately. Seal it with a zip bag and dispose of it properly.

#### 【Result interpretation】

1. When the FAM detection channel shows an S-type amplification curve, and the Ct value is < 35, the result is positive.
2. When the FAM detection channel shows no S-type amplification curve, as there is no Ct value or the Ct value is ≥35, the result is negative.

#### 【Limitations】

1. Aerosol contamination of amplification products can easily cause false positive results. The testing laboratory should be set up strictly following the requirements of the PCR testing laboratory.
2. A negative result cannot completely exclude the possibility of pathogen infection and should be combined with clinical observations and additional testing for interpretation.

#### 【Product quality indicators】

1. Sensitivity: The product detection limit of the kits is 1000 copies/ml.
2. Specificity: The test kits do not detect the cross-reaction of the pathogen samples.
3. Precision: When a strong positive sample and a weak positive sample are repeated testing for 10 times in a row separately, the CV values of their Ct values show less than 5%.

#### 【Precautions】

1. Product quality inspection: Before using PCR amplification reagent, unpack to check if the lyophilized reagent on the tube bottom is normal (white, clumped). If liquefied, it cannot be used. Otherwise, it will affect the PCR results.
2. Pipette use: When drawing 20μL of nucleic acid supernatant as the template for PCR amplification, only depress the plunger of the pipette to the first stop. Do not press to the bottom, otherwise the sample volume will exceed 20μL which will affect the results.
3. Sample information setting: Make sure that the PCR tubes are placed in the same sample well position as set in the instrument. For example, if the PCR tube is placed in the #1 well of the sample plate, then select the corresponding #1 position in the interface. After setting the sample name and the test item, press the Run icon to start the PCR amplification process.

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